

A Novel Denoising Strategy for Ultrasound Elastography



K. N. Magdoom^{1,2,3}, Julian A. Rey^{1,4} and Peter J. Basser¹


¹Eunice Kennedy Shriver National Institute of Child Health and Human Development, National Institutes of Health, Bethesda, MD, United States, ²The Military Traumatic Brain Injury Initiative (MTBI²), Uniformed Services University of the Health Sciences, Bethesda, MD, United States, ³The Henry M Jackson Foundation for the Advancement of Military Medicine (HJF) Inc., Bethesda, MD, United States, ⁴National Institute of General Medical Sciences, Bethesda, MD, United States



Disclaimer

The authors have no conflicts of interest to disclose. The views, information or content, and conclusions presented do not necessarily represent the official position or policy of, nor should any official endorsement be inferred on the part of the Uniformed Services University, the Department of Defense, the U.S. Government or the Henry M. Jackson Foundation for the Advancement of Military Medicine, Inc.

Palpation Exam



<https://en.wikipedia.org/wiki/Palpation>

Mariappan et al. (2010) Clinical Anatomy 23:497-511

US Elastography – General Principle



Palpate the sample



Measure the resulting displacement




Infer stiffness from the measured strain



$$\epsilon_l = \frac{\partial u}{\partial x}$$
$$\epsilon_{yy} = \frac{\partial u}{\partial y}$$
$$\epsilon_{eq} = \sqrt{\frac{2\epsilon_{ij}\epsilon_{ij}}{3}}$$

Noise amplification in the strain tensor calculation

Compatibility Conditions



Strain tensor

$$\frac{\partial^2 \varepsilon_{xx}}{\partial y^2} + \frac{\partial^2 \varepsilon_{yy}}{\partial x^2} = 2 \frac{\partial^2 \varepsilon_{xy}}{\partial x \partial y}$$

Enforcing the above condition on the strain tensor ensures there are **no cracks, slips, folds in the material** which the noise would mimic.

Denoising Strategy


$$\min_u \|\vec{u} - \vec{u}_{measured}\|_2$$

s.t.

$$\left| \frac{\partial^2 \varepsilon_{xx}}{\partial y^2} + \frac{\partial^2 \varepsilon_{yy}}{\partial x^2} - 2 \frac{\partial^2 \varepsilon_{xy}}{\partial x \partial y} \right| \approx 0$$

$$\varepsilon_{xx} = \frac{\partial u_x}{\partial x}, \varepsilon_{yy} = \frac{\partial u_y}{\partial y}, \varepsilon_{xy} = \frac{1}{2} \left(\frac{\partial u_x}{\partial y} + \frac{\partial u_y}{\partial x} \right)$$


Non-linear differentiation



*Essentially non-oscillatory
(ENO) scheme*

$$\frac{\partial f}{\partial x} = \begin{cases} \frac{f_{i+1} - f_i}{\Delta x} & |f_{i+1} - f_i| < |f_i - f_{i-1}| \\ \frac{f_i - f_{i-1}}{\Delta x} & \text{Otherwise} \end{cases}$$

Finite Element Simulation



Background

Young's Modulus 18 kPa
Poisson's ratio 0.49


Orange Inclusion (Stiff)

Young's Modulus 36 kPa
Poisson's ratio 0.49

Green Inclusion (Soft)


Young's Modulus 6 kPa
Poisson's ratio 0.49

Displacement Profile & Strain




Denoising Efficacy

- Additive Gaussian noise was added to the simulated displacement field to test efficacy.
- DNR – displacement noise ratio defined as ratio of maximum displacement (5 mm) and noise standard deviation
- Contrast significantly improved with denoising




Strain Elastography Measurement



- **Sample:** CIRS elasticity QA phantom (Model 049A)
- **US system :** Verasonics Vantage 128
- **Probe:** C5-2v (4 MHz center frequency)
- Pre- and post-compression B-mode images acquired using ray scanning.
- 2D displacement measured with speckle tracking using normalized cross-correlation of images.

Results



Conclusion

- We have demonstrated a novel strategy to denoise deformation fields by enforcing compatibility conditions using simulation and experimental data.
- The method is general and does not make assumptions about the material properties.
- It can be easily applied to various elastography techniques not restricted to strain elastography

Acknowledgements

This study was supported by the Intramural Research Program of the NICHD. This work was funded by the Department of Defense in the Military Traumatic Brain Injury Initiative under award, HU0001-22-2-0058.